

THE DENTAL ATTACHMENTS COMPANY







| FIXED PROSTHESIS OT BRIDGE | |
|---|------------|
| OT Equator abutment, Cuff height measurer, KIT 42, Titanium abutments | Pag. 3 |
| Castable abutments, Titanium abutments with lateral holes, Titanium abutment for endoral welding | Pag. 4 |
| Healing abutments, Impression copings | Pag. 4 |
| Laboratory analogs, Scan Abutment, Screws, OT Equator screwing tools | Pag. 5 |
| Screwing tools, Universal handle, Seeger inserter, T-Bar Components and accessories | Pag. 6 |
| Spare parts, Gingival connectors, Seeger seat cleaning cutter, Ot Cem | Pag. 7 |
| REMOVABLE PROSTHESIS ——————————————————————————————————— | |
| OT Equator for implants, OT Equator for bars, Ot Equator caps kit, OT Equator castable | Pag.8 |
| Spare parts OT Equator | Pag. 9 |
| OT Equator transfer and analogs, OT Equator tools and accessories | Pag. 10 |
| Single threaded spheres, Spare parts for single threaded spheres | Pag. 11 |
| Sphero block and Sphero flex pfor implants, Cuff height measurer | Pag. 12 |
| Tools Sphero flex and Sphero block, Anti-unscrewing system | Pag. 13 |
| Laboratory set | Pag. 14 |
| Basic laboratory set | Pag. 15 |
| OT Cap/OT Box attachments, OT Box large attachments, OT Strategy attachments, OT Bar attachments | Pag. 16 |
| OT Vertical attachments, OT Unilateral attachments, OT Cap Tecno attachments, OT Lock attachments | Pag. 17 |
| Pivots set for clinic, Ot Equator pivots set | Pag. 18 |
| Pivot OT Equator caps set, Pivot Block / Flex | Pag. 19 |
| Spare parts OT Cap normal | Pag. 20 |
| Spare parts OT Cap micro | Pag. 21 |
| OT Cap assortment sets | Pag. 22 |
| Tools and accessories OT Cap normal/micro | Pag. 23 |
| OT Strategy spare parts | Pag. 24-25 |
| OT Bar multiuse spare parts | Pag. 25 |
| OT Vertical spare parts, OT Reverse spare parts | Pag. 26 |
| OT Lock spare parts, Sphere repair sets | Pag. 27 |
| OT Lever, Small parallelometer, Denture boxes | Pag. 28 |
| OT Cem, Broken screw extractor kit | Pag. 29 |
| Rhein83: The training | Pag. 30 |
| Rhein83 direct services | Pag. 31 |

OT EOUATOR ABUTMENT



Ref. 030

• 1 OT EQUATOR TITANIUM ABUTMENT WITH TIN COATING

CUFF HEIGHT MEASURER



Ref. 008MBG

SET CONTENT

- CUFF HEIGHT SLIDER GAUGE CUFF HEIGHT FIXED ROD GAUGE
- SILICON RING DISPENSER
- SILICON RINGS

SET 42



Ref. BE42 + n° 42 OT Equator (Ref. 030 to be chosen) ASSORTMENT SET (autoclavable) + 42 OT EQUATOR +TIN ABUTMENTS THE SET42 CAN BE CUSTOMIZED ON THE BASIS OF THE IMPLANTS BRANDS, IMPLANTS DIAMETERS AND CUFF HEIGHTS DESIRED. PLEASE INDICATE, IN THE ORDER, THE QUANTITY

AND THE ABOVE SPECIFICATIONS FOR EVERY KIND OF ABUTMENT. (Tools not included in the set)



Ref. BE42

• 1 EMPTY BOX FOR ASSORTED OT EQUATOR ABUTMENTS

ABUTMENTS FOR FIXED PROSTHESIS



Ref. MTSEG

Ø 3.5 mm

Ref. MT40SEG

Ø4mm

Ref. MT50SEG

Ø5mm

EXTRAGRADE TITANIUM ABUTMENT

SET CONTENT

- EXTRAGRADE TITANIUM ABUTMENT TO BE CEMENTED h9,5mm TITANIUM SCREW h2mm
- WHITE SEEGERS WITH HANDLE PINK SEEGER WITH HANDLE



Ref. MMSEG

Ø3.5 mm

Ref. MM40SEG

MONCONE IN TITANIO EXTRAGRADE MINI

SET CONTENT

- MINI EXTRAGRADE TITANIUM ABUTMENT TO BE CEMENTED h5,7mm TITANIUM SCREW FOR MINI EXTRAGRADE TITANIUM ABUTMENT
- WHITE SEEGERS WITH HANDLE PINK SEEGER WITH HANDLE



Ref. MTSFEG

Ø 3.5 mm

Ref. MTEG15SF

EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15°

SET CONTENT

- EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15° TO BE CEMENTED
- WHITE SEEGERS WITH HANDLE
- PINK SEEGER WITH HANDLE SPHERICAL TITANIUM SCREW FOR 15° INCLINED ABUTMENT

ABUTMENTS FOR FIXED PROSTHESIS



Ref. CMEG40C

EXTRAGRADE CASTABLE ABUTMENT

SET CONTENT

- EXTRAGRADE CASTABLE ABUTMENT WITH HOLE TITANIUM SCREW h2 mm
- WHITE SEEGERS WITH HANDLE PINK SEEGER WITH HANDLE



Ref. MCEG15SF

EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE

SET CONTENT

- EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE (ONLY AVAILABLE WITH Ø4mm)
- 2 WHITE SEEGERS WITH HANDLE
- PINK SEEGER WITH HANDLE



Ref. CMTFP

h 1.8 mm

Ref. CMTB30FP

h3mm

TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCMENT

SET CONTENT

- TITANIUM ABUTMENT WITH LATERAL HOLES (ONLY AVAILABLE WITH Ø 3,5mm)
- TITANIUM SCREW h2 mm WHITE SEEGERS WITH HANDLE
- PINK SEEGER WITH HANDLE



Ref. MT40GL

TITANIUM ABUTMENT FOR ENDORAL WELDING

SET CONTENT

- TITANIUM NORMAL ABUTMENT FOR INTRA ORAL WELDING h9,5 mm Ø4mm TITANIUM SCREW h2 mm
- WHITE SEEGERS WITH HANDLE
- PINK SEEGER WITH HANDLE

HEALING ABUTMENTS



Ref. 159VAG

SCREW AND HEALING RING

SET CONTENT

- TITANIUM SCREW
- TITANIUM HEALING RING



Ref. 159EAG

OT EQUATOR AND HEALING RING

- SET CONTENT
 1 TITANIUM + TIN OT Equator SCREW FOR HEALING RING
- TITANIUM HEALING RING



Ref. 159VEA

Titanium + TiN OT Equator screw (for self-extracting peek ring or titanium self-extracting healing ring)

IMPRESSION COPINGS



Ref. 144TLE

• 1 TITANIUM OT Equator LONG IMPRESSION COPING + TITANIUM LONG SCREW h11,5mm



Ref. 144TTE

• 1 TITANIUM OT Equator IMPRESSION COPING + TITANIUM SCREW h6 mm



Ref. 144MTE

• 2 OT Equator IMPRESSION COPINGS FOR INDIVIDUAL TRAY



Ref. TPM2

• 2 MINI OT Equator IMPRESSION COPINGS h3,7mm

LABORATORY ANALOGS



Ref. 144AE

• 2 LABORATORY ANALOGS ø 4mm

Ref. AI502

• 2 LABORATORY ANALOGS ø 5mm



Ref. 144AVC4

- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 4mm
- 1 LABORATORY ANALOG WITH SCREW FOR CAD/CAM ø 5mm

OT EQUATOR SCAN ABUTMENT



Ref. 145SAE

• 1 OT Equator TITANIUM SCAN ABUTMENT + TITANIUM SCREW

SCREWS



Ref. 147VME

• 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT (for screwdriver/stem for inclined Channel))

Ref. 146VMI

• 1 SPHERICAL TITANIUM SCREW FOR INCLINED ABUTMENT (for screwdriver/stem for inclined Channel))

•

Ref. 148VCM

 TITANIUM SCREW FOR MINI EXTRAGRADE ABUTMENT (for screwdriver/stem for inclined Channel))

Ref. VC

• 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h2mm



Ref. VM

- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h6mm
- Ref. VXL
- 1 TITANIUM SCREW FOR EXTRAGRADE ABUTMENT h11,5mm

TOOLS



Ref. 774CHE

Ref. 775CLE*

(*Long version)

 SQUARE SCREW DRIVER + HOLDER for tightening OT Equator (SQUARE TIP 1.25MM) (USABLE WITH THE MANUAL TORQUE WRENCH)



Ref. 774HC

• 1 INTERCHANGEABLE HOLDER



Ref. 760CSC


Ref. 760CS*

(*Long version)

1 SCREWDRIVER/STEM FOR ANGULATED CHANNEL FOR TORQUE CONTROLLER CONTRA ANGLE (SQUARE 1.25MM) AND ADAPTER FOR MANUAL TORQUE WRENCH

FIXED PROSTHESIS OT BRIDGE

TOOLS •



MEDICAL STEEL T-BAR

Ref.190IS



 4 MEDICAL STEEL JOINTS Ref. G80A

 4 MEDICAL STEEL CYLINDERS Ref. T20A

• 4 MEDICAL STEEL RETENTIVE CYLINDERS Ref. T20Ai

TITANIUM T-BAR •



Ref. G80T 4 TITANIUM JOINTS

 4 TITANIUM CYLINDERS Ref. T20T

 4 TITANIUM RETENTIVE CYLINDERS Ref. T20Ti

TOOLS T-BAR



Ref. T-APLUS 2 T-BAR ADAPTERS

SPARE PARTS —

| | Ref. RMTSFEG | • 1 EXTRAGRADE TITANIUM ABUTMENT INCLINED AT 15° Ø 3,5mm | |
|--|----------------|---|--|
| | Ref. RMTEG15SF | • 1 EXTRAGRADE TITANIUM ABUTMENT INCLINED AT15° Ø 4mm | |
| <u>a</u> | Ref. RMTSEG | • 1 EXTRAGRADE TITANIUM ABUTMENT Ø 3,5mm | |
| TE. | Ref. RMT40SEG | • 1 EXTRAGRADE TITANIUM ABUTMENT Ø 4mm | |
| | Ref. RMT50SEG | • 1 EXTRAGRADE TITANIUM ABUTMENTØ 5mm | |
| | Ref. RMMSEG | • 1 MINI EXTRAGRADE TITANIUM ABUTMENT Ø 3mm | |
| (3) | Ref. RMM40SEG | • 1 MINI EXTRAGRADE TITANIUM ABUTMENT Ø4mm | |
| | D. C. COMP. | 1. CASTADLE SUFFEYE FOR INCLINED TITANIUM ARRITMENT WITHOUT HOLE GOS | |
| B.A. | Ref. GCMR | • 1 CASTABLE SLEEVE FOR INCLINED TITANIUM ABUTMENT WITHOUT HOLE Ø 3,5mm | |
| 1 | Ref. GCMR4 | • 1 CASTABLE SLEEVE FOR INCLINED TITANIUM ABUTMENT WITHOUT HOLE Ø4mm | |
| | Ref. GC | • 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø 3,5mm | |
| 0 | Ref. GC4 | • 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø4mm | |
| | Ref. GC5 | • 1 CASTABLE SLEEVE FOR TITANIUM ABUTMENT Ø 5mm | |
| | Ref. GMM | • 1 CASTABLE SLEEVE FOR MINI TITANIUM ABUTMENT Ø 3,5mm | |
| 0 | Ref. GMM4 | • 1 CASTABLE SLEEVE FOR MINI TITANIUM ABUTMENT Ø 4mm | |
| Eb. | Ref. MTFP | • 1 TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE WITH WIRE-REINFORCMENT h1,8 mm (ONLY AVAILABLE WITH Ø 3,5 mm) | |
| | Ref. MTB30FP | TITANIUM ABUTMENT WITH LATERAL HOLES FOR TEMPORARY RESIN BRIDGE | |
| | nei ini baari | WITH WIRE-REINFORCMENT h3 mm (ONLY AVAILABLE WITH Ø 3,5mm) | |
| | Ref. HR800 | • 1 TITANIUM WIRE WITHOUT ELASTIC MEMORY Ø 1,0mm | |
| | Ref. MEG40C | • 1 EXTRAGRADE CASTABLE ABUTMENT Ø 4mm | |
| 9 | Ref. RMCEG15SF | • 1 EXTRAGRADE CASTABLE ABUTMENT WITHOUT HOLE (ONLY AVAILABLE WITH Ø 4,0mm)) | |
| | Ref. SM6 | • 6 WHITE SEEGERS WITH HANDLE (STANDARD RETENTION) | |
| 7 | Ref. SRM6 | • 6 PINK SEEGER WITH HANDLE (SOFT RETENTION) | |
| | Ref. 022RGO | • 3 CASTABLE CONNECTION BAR | |
| | Ref. FSS | • 1 BUR FOR SEEGER HOUSING POLISHING | |
| The state of the s | Ref. OC | • 1 OT CEM (COMPOSITE CEMENT: 1 Two-component syringe , 10 AUTOMIX TIPS) | |
| | Ref. TOC | • 10 AUTOMIX TIPS | |

Ref. 130

OT EQUATOR FOR IMPLANTS =



OT EQUATOR FOR IMPLANTS SET

SET CONTENT

- OT EQUATOR TITANIUM + TIN ABUTMENT
- STAINLESS STEEL HOUSING PROTECTIVE DISK
- RETENTIVE CAPS
 - (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



OT EQUATOR AND SMARTBOX FOR IMPLANTS SET

SET CONTENT

- OT EQUATOR TITANIUM + TIN ABUTMENT
- ANODIZED SMART BOX HOUSING WITH BLACK POSITIONING CAP
- PINK PROTECTIVE DISK
 - (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)



Ref. 030

Ref. 131

OT EQUATOR TITANIUM + TIN ABUTMENT

THE OT EQUATOR ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING INDICATE: IMPLANT, DIAMETER AND CUFF HEIGHT. CUFF AVAILABLE TO 7mm.

OT EQUATOR FOR BARS -



OT EQUATOR THREADED SLEEVE SYSTEM SET FOR BONDING

SET CONTENT

- EQUATOR THREADED for titanium sleeve (1.6 mm THREAD)
- THREADED SLEEVES (1.6 mm THREAD) STAINLESS STEEL HOUSINGS • 2
- OT EQUATOR SPACERS RETENTIVE CAPS

(2 EXTRA-SOFT, 2 SOFT, 2 STANDARD, 2 PROCESSING)

ELASTIC SEEGER PASSIVE BAR CONNECTIONS SYSTEM



Ref. 158ESA

Ref. 160EQB

ELASTIC SEEGER SET

- **SET CONTENT**
- CASTABLE CYLINDER FOR SEEGER
- SELF-EXTRACTING SEEGER
 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER

CAPS ASSORTMENT SET



Ref. 192ECE

CAPS ASSORTMENT SET SET CONTENT

- STAINLESS STEEL HOUSING
- RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)
- PROTECTIVE DISK

SMART BOX SET SELF-ALIGNING OT EQUATOR HOUSING



Ref. 335SBC

SMART BOX SET



- SET CONTENT SMART BOX HOUSING WITH BLACK POSITIONING CAP
- PINK PROTECTIVE DISK
- RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)

OT EQUATOR CASTABLE —



Ref. 092ECQ

OT EQUATOR CASTABLE SET

SET CONTENT

- SINGLE MALES STAINLESS STEEL HOUSINGS
- RETENTIVE CAPS (2 SOFT, 2 STANDARD)

SPARE PARTS OT EQUATOR —

| | | - | | | |
|-----|---|----------------|-----|--|--|
| 9 | | Ref. 140CEV | • 4 | RETENTIVE CAPS - VIOLET "STRONG" 2700g | |
| Q | | Ref. 140CET | • 4 | RETENTIVE CAPS - CLEAR "STANDARD" 1800g | |
| | | Ref. 140CER | • 4 | RETENTIVE CAPS - PINK "SOFT" 1200g | |
| | | Ref. 140CEG | • 4 | RETENTIVE CAPS - YELLOW (EXTRA-SOFT 600g) | |
| 9 | | Ref. 140CEN | • 4 | PROCESSING CAPS - BLACK (FOR LABORATORY USE ONLY) | |
| 4 | | Ref. 140CEDR8 | • 4 | RETENTIVE CAPS - AQUA CAPS UNDERSIZED INTERNAL DIAMETER 2,2mm (STRONG 1300G) | |
| 4 | | Ref. 140CEADR8 | • 4 | RETENTIVE CAPS - ORANGE UNDERSIZED INTERNAL DIAMETER 2,2mm (VERY ELASTIC 350G) | |
| 1 | | Ref. 141CAE | • 2 | STAINLESS STEEL HOUSINGS | |
| 1 | | Ref. 141CTE | • 2 | TITANIUM HOUSINGS FUCHSIA ANODIZED | |
| | | Ref. 330SBE | • 1 | SMART BOX HOUSING (WITH BLACK POSITIONING CAP) | |
| | | Ref. 335CSB | • 4 | BLACK SMARTBOX (POSITIONING CAPS) | |
| 4 | | Ref. 039SFE2 | • 1 | OT EQUATOR THREADED 2mm (DIRECT INTO 2mm THREAD) | |
| - 4 | | Ref. 339SFE | • 1 | EQUATOR THREADED 1.6mm (FOR TITANIUM SLEEVE) | |
| | 3 | Ref. 239GSF | • 1 | TITANIUM SLEEVE (1.6mm THREAD) | |
| 1 | | Ref. 339SSE | •1 | OT EQUATOR SPACER | |
| 4 | 3 | Ref. 159VEA | •1 | OT EQUATOR SEEGER (FOR SELF-EXTRACTING PEEK RING OR TITANIUM HEALING RING) | |
| * | | Ref. 158VAE | • 1 | TITANIUM SCREW FOR SELF-EXTRACTING SEEGER | |
| | n | Ref. 158SAE | • 1 | SELF-EXTRACTING WHITE SEEGER (FOR BAR LOCKING) | |
| | | Ref. 158CCE | •1 | CASTABLE CYLINDER height 2,5mm (FOR SEEGER) | |
| | | Ref. 158CCAE | •1 | HIGH CASTABLE CYLINDER height 3,5mm (FOR SEEGER) | |
| 4 | 3 | Ref. 151SS | • 2 | CASTABLE MALES OT EQUATOR | |
| | | | | | |

ACCESSORIES OT EQUATOR



Ref. 144TLE • 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW




Ref. 144TTE • 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW



Ref. 144MTE • 2 IMPRESSION COPINGS (FOR INDIVIDUAL TRAY)



IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR Ref. 044CAIN





Ref. 144AE • 2 LABORATORY ANALOGS ø 4mm

• 2 LABORATORY ANALOGS ø 5mm **Ref. AI502**

• 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø4mm **Ref. 144AVC4**

Ref. 144AVC5 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø5mm



Ref. 161SCME OT EQUATOR SPACERS FOR HOUSING

EQUATOR TOOLS



Ref. 74AC01 PARALLELOMETER KEY "NORMAL SIZE"



OT FOUATOR SOLIARE SCREW DRIVER + HOLDER Ref. 774CHE

(SQUARE 1.25MM) (USABLE WITH RATCHET TORQUE CONTROL DEVICE)



OT EQUATOR INTERCHANGEABLE HOLDER



SQUARE DRIVER CONNECTOR FOR TORQUE **Ref. 760CE** CONTROLLER CONTRA ANGLE (SQUARE 1.25mm)

Ref. 485IC

MULTIUSE CAPS INSERTER/EXTRACTOR TOOL • 1



Ref. 185SIS METAL INSERTION TOOL FOR SEEGER USE TOGETHER WITH 491EC

MULTIUSE EXTRACTOR TOOL FOR CAPS **Ref. 491EC** (WITH HOLE STEM HOLDER)

RATCHET TOROUE CONTROL DEVICE (FOR OT EOUATOR, SPHERO BLOCK Ref. 760CRD-US

SINGLE THREADED SPHERES WITH THREADED SLEEVE TO BE BONDED NORMAL - MICRO



SET SINGLE THREADED SPHERE NORMAL SIZE + TIN WITH SLEEVE TO BE BONDED SET CONTENT

- TITANIUM NORMAL SINGLE THREADED SPHERE • 2 (HEX 1,3mm THREAD 1,6mm)
- TITANIUM BOND-IN THREADED SI FEVE
- S/STEEL SPACER for the THREADED SLEEVE (NORMAL sphere)



SET SINGLE THREADED SPHERE MICRO SIZE + TIN WITH SLEEVE TO BE BONDED

SET CONTENT

- TITANIUM MICRO SINGLE THREADED SPHERE (HEX 0,9mm THREAD 1,6mm)
- TITANIUM BOND-IN THREADED SLEEVE
- S/STEEL SPACER for the THREADED SLEEVE (NORMAL sphere)

SPARE PARTS



Ref. 239SFN

Ref. 139KSFN

Ref. 139KSFM

SINGLE THREADED SPHERE NORMAL SIZE + TiN FOR THREADED SLEEVE HEX 1,3mm THREAD 1,6mm



Ref. 239SFM



Ref. 239GSF

SINGLE THREADED SPHERE MICRO SIZE + TIN FOR THREADED SLEEVE HEX 0,9mm THREAD 1,6mm



TITANIUM BOND-IN THREADED SLEEVE FOR NORMAL AND MICRO - 1,6mm



Ref. 239SSN

• 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (NORMAL sphere)



Ref. 239SSM

• 1 STAINLESS STEEL SPACER FOR THREADED SLEEVE (Micro sphere)



Ref. 039SFN2

SINGLE THREADED SPHERE NORMAL SIZE + TIN HEX 1.3mm THREAD 2mm



Ref. 039SFM2

SINGLE THREADED SPHERE MICRO SIZE + TIN HEX 0,9mm THREAD 2mm



Ref. 760CBR

HEX DRIVER FOR THREADED SPHERE NORMAL HEX 1,27mm



Ref. 772CSF

HEX DRIVER FOR THREADED SPHERE MICRO HEX 0,9mm

OT CAP IMPRESSION COPINGS



Ref. 044CAIN

IMPRESSION COPINGS NORMAL SIZE Ø2,5mm WITH INTERCHANGEABLE CAP



Ref. 044CAI22

IMPRESSION COPINGS • 2 Ø2,25mm SPHERES WITH INTERCHANGEABLE CAP



Ref. 044CAIM

• 2 IMPRESSION COPINGS MICRO SIZE Ø1,8mmWITH INTERCHANGEABLE CAP

IMPLANTOLOGY TITANIUM ATTACHMENTS FOR OVERDENTURE —

SPHERO FLEX - SPHERO BLOCK

Ref. 109



SPHERO FLEX TIN COATED

SET CONTENT

- TITANIUM ABUTMENT WITH SWIVEL SPHERE
- (SPHERE Ø 2,5mm) CAPS (Yellow Pink Clear) STAINLESS STEEL HOUSING
- 1
- PROTECTIVE DISK DIRECTIONAL RINGS

Ref. 002



SPHERO BLOCK NORMAL TIN COATED

SET CONTENT

- TITANIUM ABUTMENT
- (SPHERE Ø 2,5 mm) CAPS (Yellow Pink Clear)
- STAINLESS STEEL HOUSING PROTECTIVE DISK
- DIRECTIONAL RINGS
- **Ref. 003**



SPHERO BLOCK MICRO TIN COATED

SET CONTENT

- TITANIUM ABUTMENT (SPHERE Ø 1,8mm)
- CAPS (Yellow Pink Clear) STAINLESS STEEL HOUSING PROTECTIVE DISK • 3
- DIRECTIONAL RINGS

THE SPHERO FLEX, SPHERO BLOCK ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS.UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN (B.D.G) mm ORDERING, PLEASE ALWAYS INDICATE: IMPLANT BRAND, IMPLANT DIAMETER AND CUFF HEIGHT. FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 0,5 TO 7mm.



CUFF HEIGHT MEASURER


THIS TOOL IS USED TO DETERMINE THE **CUFF HEIGHT OF IMPLANTS WITH** INTERNAL OR EXTERNAL HEXAGON

Ref. 008MBG




SET CONTENT

- CUFF HEIGHT SLIDER GAUGE
- CUFF HEIGHT FIXED ROD GAUGE SILICON RING DISPENSER
- SILICON RINGS



TOOLS SPHERO FLEX - SPHERO BLOCK



ELASTIC DOWEL

ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUEST TO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

Extra charge for the dowel application



ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR OT VERTICAL - OT TECNO - TITANIUM SINGLE SPHERE

MASTER STARTER SET FOR LABORATORY



SET CONTENT

OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
- castable spheres (NORMAL and micro size) titanium single spheres (NORMAL and micro size)
- + sphere holder + 2 sliding base castable OT CAP bars (NORMAL and micro size) complete castable OT CAP TECNO bars (NORMAL and micro size)
- castable OT BOX CLASSIC bars (upper + lower part)
- NORMAL and micro size
- castable OT BOX SPECIAL bars (NORMAL and micro size) castable OT BOX connectors castable OT BOX MONO housings (NORMAL and micro size)
- plastic positioner rings (NORMAL and micro size)
- OT CAP retentive caps (NORMAL and micro size) Clear- standard retention - PINK soft retention YELLOW extra soft retention - GREEN very elastic retention
- TITAN CAP titanium caps (NORMAL and micro size)
 BLACK CAPS- only for lab use (NORMAL and micro size)
- S/STEEL housings • 8
- ANALOGUES for plaster models (NORMAL and micro size)
- protective disks
- aluminium spacers (NORMAL and micro size)

OT STRATEGY:

- castable spheres (2 standard and 2 high)
- castable STEADY
- retentive caps: 2 clear 2 pink 2 yellow for housing
- black caps for housing only for lab use
- s/steel housings
- plastic positioner rings
- retentive caps for duplication technique: CLEAR - standard retention - PINK soft retention YELLOW extra soft retention
- black caps for duplication technique only for lab use
- ANALOGUES for plaster models

OT BAR MULTIUSE:

- castable bar • 2
- gingival connector
- positioner clip A
- positioner clip B
- castable housings
- clips 2 PINK soft retention 2 YELLOW medium retention

Ref. 005SKLUS

OT VERTICAL:

- castable male • 2
- castable STEADY
- 4 clips - 2 WHITE- standard rention - 2 GREEN soft retention
- castable parallelometer keys
- ceramic bars

TOOLS:

- 2
- Pincers 2 Reamer tools for caps (NORMAL and micro size)
 Reamer tools for OT STRATEGY caps
 Parallelometer keys for OT CAP (NORMAL and micro size)
 Parallelometer key for OT STRATEGY size
 Parallelometer keys for OT CAP TECNO (NORMAL and micro size)
 Parallelometer keys for OT BAR MULTIUSE
 Blue Plastic "multiuse" insertion tool
- (for steel tool and for inserting the caps in the mouth)
- OT CAP insertion tool
- OT STRATEGY insertion tool
- OT BAR clips insertion tool
- OT VERTICAL clips insertion tool

ATTACHMENTS OT CAP & OT BOX - OT STRATEGY- OT BAR OT VERTICAL - OT UNILATERAL - OT EQUATOR

BASIC LABORATORY SET



OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (NORMAL and micro size)
- castable spheres NORMAL/micro
- castable OT CAP bars (NORMAL and micro size) castable OT BOX CLASSIC bars (upper + lower part) NORMAL and micro size castable OT BOX SPECIAL bars (NORMAL and micro size)
- castable OT BOX connectors
- castable OT BOX MONO housings (NORMAL and micro size)
- plastic positioner rings (NORMAL and micro size)
- OT CAP assorted retentive caps (NORMAL and micro size) PINK soft retention - YELLOW extra soft retention
- 13 BLACK CAPS- only for lab use (NORMAL and micro size)
 4 S/STEEL housings (NORMAL and micro size)

- castable spheres (2 standard and 2 high) castable STEADY
- pink retentive caps
- s/steel housings
- plastic positioner rings
- retentive caps for duplication technique: PINK soft retention YELLOW extra soft retention
- 4 black caps for duplication technique only for lab use

OT EQUATOR:

- single males
- stainless steel housings
- retentive caps (2 soft, 2 standard)
- processing black caps

OT UNILATERAL:

- unilateral castable attachment
- uni-box castable
- plastic positioner micro size
- ot cap micro size: 1 pink 1 black processing cap
- ot strategy (for duplication technique): 1 yellow 1 black processing cap
- connector

OT BAR MULTIUSE:

- castable bar
- gengival connector positioner clip (A - B)

- clips 2 PINK soft retention 2 YELLOW medium retention

OT VERTICAL:

Ref. 005SKLBUS

- castable male
- castable STEADY
- clips 2 WHITE- standard rention 2 GREEN soft retention
- castable parallelometer keys
- ceramic bars

TOOLS:

- Tweezer
- Parallelometer keys for OT CAP NORMAL
- Parallelometer keys for OT CAP micro
- Parallelometer key for OT STRATEGY
- Parallelometer key for OT BAR MULTIUSE
- . 1 Blue Plastic "multiuse" insertion tool (for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- OT STRATEGY insertion tool • 1
- OT BAR clips insertion tool
- OT VERTICAL clips insertion tool

OT CAP AND OT BOX ATTACHMENTS •



NORMAL

Ref. 092CAN

SET CONTENT

OT CAP NORMAL/MICRO CASTABLE SINGLE SPHERES

MICRO

PINK RETENTIVE CAPS STAINLESS STEEL HOUSINGS PLASTIC POSITIONERS

Ref. 092CAM

NORMAL

MICRO

Ref. 196BCN

OT CAP NORMAL/MICRO "ECONOMIC"

SET CONTENT

CASTABLE PRESECTIONED BAR • 1

CASTABLE PRESECTIONED BEVELED BAR
CLEAR RETENTIVE CAPS (Standard retention)

NORMAL

Ref. 099BSN

Ref 197BCM

OT CAP & OT BOX MONO NORMAL/MICRO SET CONTENT

CASTABLE PRESECTIONED BARS (1 straight + 1 beveled) CLEAR RETENTIVE CAPS

CASTABLE MONO OT BOX

PLASTIC POSITIONERS

NORMAL

Ref. 058BSN

Ref. 058BSM

Ref. 153BCN

Ref. 099BSM

SET CONTENT

MICRO

OT BOX SPECIAL BARS PLASTIC POSITIONERS • 2 • 4

CONNECTORS

NORMAL

MICRO

OT BOX CLASSIC NORMAL/MICRO + CONNECTORS

OT BOX SPECIAL NORMAL/MICRO + CONNECTORS

SET CONTENT

• 2 LOWER BARS

PLASTIC POSITIONERS

Ref. 153BCM

• 4 CONNECTORS

OT BOX LARGE OVERSIZED CASTABLE HOUSING FOR CHAIRSIDE POSITIONING OF CAPS

NORMAL Ref. 253BLN

OT BOX LARGE NORMAL/MICRO + CONNECTORS SET CONTENT

MICRO

• 2

Ref. 253BLM

OT BOX LARGE BARS PLASTIC POSITIONERS

CONNECTORS

OT STRATEGY ATTACHMENTS -



Ref. 098SSS

OT STRATEGY CAPS FOR DUPLICATION TECHNIQUE

SET CONTENT

CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH) CASTABLE STEADY

• 2

RETENTIVE CAPS FOR DUPLICATION TECHNIQUE



Ref. 098CAL

OT STRATEGY CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE SET CONTENT

CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)

RETENTIVE CAPS FOR HOUSING TECHNIQUE (2 SOFT, 2 STANDARD) STAINLESS STEEL HOUSINGS

PLASTIC POSITIONERS

CASTABLE STEADY

OT BAR MULTIUSE ATTACHMENTS =



Ref. 0210BM

OT BAR MULTIUSE + CONNECTORS

SET CONTENT
• 2 CASTABLE BARS

POSITIONER CLIP A
POSITIONER CLIP B

. 4 CASTABLE BOXES

RETENTIVE PINK CLIPS • 4

RETENTIVE YELLOW CLIPS RETENTIVE RED CLIPS

CONNECTORS GINGIVAL CONNECTOR

OT VERTICAL ATTACHMENTS



Ref. 0710BV

OT VERTICAL

SET CONTENT

- CASTABLE MALES • 4
- 2
- CASTABLE STEADY
 RETENTIVE WHITE CLIPS
 RETENTIVE GREEN CLIPS
- CASTABLE PARALLELOMETER KEYS + PIN

OT UNILATERAL EXTRA CORONAL ATTACHMENTS WITH COMBINED FUNCTION

Ref. 098UNL

OT UNILATERAL + CONNECTOR

SET CONTENT

- UNILATERAL CASTABLE ATTACHMENT UNI-BOX CASTABLE
- PLASTIC POSITIONER MICRO SIZE PINK CAP MICRO SIZE

- STRATEGY YELLOW CAP (for duplication technique) BLACK CAP MICRO SIZE (for laboratory use only)
- STRATEGY BLACK CAP (for laboratory use only)
- CONNECTOR



Ref. 156BCU

• 1 UNILATERAL CASTABLE ATTACHMENT



Ref. 156CCU

• 1 UNI-BOX CASTABLE

OT CAP TECNO

NORMAL


Ref. 093CTN


OT CAP TECNO NORMAL/MICRO SET CONTENT


PRE-ANGULATED CASTABLE EXTENSIONS TITANIUM SINGLE THREADED SPHERES

Ref. 093CTM

TITANIUM BOND-IN THREADED SLEEVES TRANSPARENT CAPS NORMAL/MICRO







Ref. 193BCA

• 2 PRE-ANGULATED CASTABLE EXTENSIONS

OT LOCK LOCKING PIN TITANIUM AND CASTABLE



Ref. 880CLT

OT LOCK SET TITANIUM

SET CONTENT

- LOCKING PIN SET
- HOUSING SHAPER CERAMIC PIN



Ref. 880CLR

OT LOCK SET TITANIUM ADJUSTABLE

SET CONTENT

- LOCKING PIN SET
- HOUSING SHAPER CERAMIC PIN
- CASTABLE SPACER RINGS

OVERDENTURE ATTACCHMENTS -

S.P.L. PIVOTS - COPING COVER

PROMOTIONAL SET FOR CLINIC



Ref. 011SPL

SET CONTENT

- 24 CASTABLE ASSORTED PIVOTS MICRO SIZE
- 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE
- 12 TITANIUM MICRO PIVOTS
- 12 TITANIUM NORMAL PIVOTS
- 16 NORMAL CAPS (Clear standard retention PINK soft retention)
- 16 MICRO CAPS (Clear standard retention PINK soft retention)
- 20 PROTECTIVE DISKS
- 2 REGULATING TOOLS (1 NORMAL and 1 micro size)
- BLUE PLASTIC "MULTIUSE" INSERTION TOOL • 1 (FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)
- 1 OT CAP insertion tool
- 4 NORMAL STAINLESS STEEL HOUSING4 MICRO STAINLESS STEEL HOUSING
- 12 DIRECTIONAL RINGS (grey 0° green 7° red 14°)
- 1 TWEEZER



Lenghts available: 10 - 9 - 7mm

Ref. 333PEK7

Ref. 333PEK9

Ref. 333PEK10

OT EQUATOR TITANIUM TIN COATED PIVOTS SET

LOW PROFILE HEAD PIVOTS

SET CONTENT

- TITANIUM OT EQUATOR PIVOTS h7,9 or 10mm (of the same length chosen)
 OT EQUATOR ANODIZED TITANIUM HOUSINGS
 YELLOW EXTRA SOFT CAPS
- 2
- PINK SOFT CAPS CLEAR STANDARD CAPS
- TRANSPARENT PROTECTIVE DISKS

Lenghts available: 10 - 9 - 7mm

Ref. 336PTE7

Ref. 336PTE9

Ref. 336PTE10

OT EQUATOR TITANIUM TIN COATED PIVOTS LOW PROFILE HEAD PIVOTS SET CONTENT

• 1 TITANIUM OT EQUATOR PIVOT h7,9 or 10mm (of the length chosen)



OT EQUATOR CAPS ASSORTMENT SET

SET CONTENT

- STAINLESS STEEL HOUSING • 1
- PROCESSING CAP RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)
- PROTECTIVE DISK



Ref. 033PSF7R

Ref. 033PSF9R

Ref. 033PSF10R

PIVOT FLEX TIN COATED

PTITANIUM PIVOT WITH SWIVEL SPHERE - NORMAL size Adapted for COPING COVER **SET CONTENT**

- TITANIUM PIVOT WITH SWIVEL SPHERE
- STAINLESS STEEL HOUSING PINK CAPS NORMAL Size Soft retention
- PROTECTIVE DISK • 1
- **DIRECTIONAL RINGS**

Lenghts available:: 10 - 9 - 7mm NORMAL ø2,5 Micro ø1,8

Ref. 036PTN7

Ref. 036PTN9

Ref. 036PTN10

TITANIUM PIVOTS SET NORMAL SIZE

SET CONTENT

TITANIUM PIVOTS (Adapted for COPING COVER) • 5



Ref. 036PTM7

Ref. 036PTM9

Ref. 036PTM10

TITANIUM PIVOTS SET MICRO SIZE

SET CONTENT

TITANIUM PIVOTS (Adapted for COPING COVER)

NORMAL

Ref. 192ACN

MICRO

CAPS ASSORTMENT SET OT CAP NORMAL/MICRO SET CONTENT



Ref. 192ACM





STAINLESS STEEL HOUSING

PROCESSING CAP RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)

PROTECTIVE DISK



Ref. 034PRK

OT REVERSE 3 OVERDENTURE POSTS FOR ROOTS, AVAILABLE IN ONE STANDARD SIZE

SET CONTENT

ROOT POSTS IN TITANIUM + TIN RETENTIVE MALE ATTACHMENT TITANIUM + NYLON

PLASTIC MANUAL TOOLS STAINLESS STEEL HOUSINGS

PINK PROTECTIVE DISKS

Ø del perno sottosfera:



Ref. 010PSP

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

SET CONTENT

20 CASTABLE PIVOTS NORMAL SIZE

Lenghts available: 14 - 12 - 10 - 9 - 7mm



Ref. 012PSM

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

SET CONTENT

• 20 CASTABLE PIVOTS MICRO SIZE

Lenghts available: 10 - 9 - 7mm



Ref. A01MOG

MOOSER BUR FOR PIVOTS OF 7 - 9 - 10mm



Ref. A03MOB

MOOSER BUR FOR PIVOTS OF 12 - 14mm



Ref. 121FD

OT REVERSE 3 DIAMOND SHAPED BUR FOR ROOTS



Ref. 100AD

DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK COLOURED RINGS WITH INCLINATION OF 0° - 7°

Ref. 485IC

• 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - NORMAL - MICRO

SPARE PARTS OT CAP SYSTEM NORMAL Ø 2,5 mm Ref. 055SCN 4 CASTABLE SINGLE SPHERES NORMAL SIZE • 2 CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE **Ref. 151BAN** • 2 CASTABLE PRESECTIONED BARS NORMAL SIZE Ref. 150BPN SILVER CAPS NORMAL SIZE (EXTRA-RESILIENT) Ref. 048CAN (ELASTIC AND GUMMY) - 350g • 6 GOLD CAPS NORMAL SIZE (EXTRA-RESILIENT) Ref. 048CON (SLIGHTLY ELASTIC) - 500g **GREEN CAPS NORMAL SIZE** Ref. 049PCN (VERY ELASTIC RETENTION) - 350g YELLOW CAPS NORMAL SIZE Ref. 060CRN AY • 6 (EXTRA SOFT RETENTION) - 500g PINK CAPS NORMAL SIZE Ref. 040CRN SN • 6 (SOFT RETENTION) - 900g Ref. 040CRN CLEAR CAPS NORMAL SIZE (STANDARD RETENTION) - 1300g • 2 TITAN CAPS NORMAL SIZE Ref. 040TCN • 6 BLACK CAPS NORMAL SIZE (FOR LABORATORY USE ONLY) Ref. 043CLN ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 350g Ref. 049PCNDR8 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 500g Ref. 060CRNAYDR8 PINK CAPS NORMAL SIZE (SOFT RETENTION) Ref.040CRNSNDR8 (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE) - 900g AQUA CAPS NORMAL SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø 2.25 SPHERE)) - 1300g Ref. 040CRNDR8 MONO OT BOX NORMAL SIZE Ref. 042AMN (CASTABLE HOUSINGS FOR CAPS) Ref. 041CAN 2 STAINLESS STEEL HOUSINGS NORMAL SIZE 2 TITANIUM HOUSINGS NORMAL SIZE Ref. 041CTN • 2 PLASTIC POSITIONERS NORMAL SIZE Ref. 101PPN NORMAL SPHERE ANALOGS Ref. 044PPN • 2

Ref. 100SAN

• 10 MINIUM SPACERS FOR OT BOX BARS NORMAL SIZE

SPARE PARTS OT CAP SYSTEM MICRO Ø 1,8 mm

| | Ref. 055SCM | • 4 | CASTABLE SPHERES MICRO SIZE | |
|-----------|----------------|------|--|--|
| Sand Sand | Ref. 151BAM | • 2 | CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE | |
| No. | Ref. 150BPM | • 2 | CASTABLE PRESECTIONED BARS MICRO SIZE | |
| | Ref. 048CAM | • 6 | SILVER CAPS MICRO SIZE (EXTRA-RESILIENT) (ELASTIC AND GUMMY) - 200g | |
| | Ref. 048COM | • 6 | GOLD CAPS MICRO SIZE (EXTRA-RESILIENT) (SLIGHTLY ELASTIC) - 450g | |
| | Ref. 049PCM | •6 | GREEN CAPS MICRO SIZE (VERY ELASTIC RETENTION) - 200g | |
| | Ref. 060CRM AY | •6 | YELLOW CAPS MICRO SIZE (EXTRA SOFT RETENTION) - 450g | |
| | Ref. 040CRM SN | •6 | PINK CAPS MICRO SIZE (SOFT RETENTION) - 800g | |
| | Ref. 040CRM | • 6 | CLEAR CAPS MICRO SIZE (STANDARD RETENTION) - 1100g | |
| | Ref. 040TCM | • 2 | TITAN CAPS MICRO SIZE | |
| | Ref. 043CLM | •6 | BLACK CAPS MICRO SIZE (FOR LABORATORY USE ONLY) | |
| | Ref. 049PCMDR8 | • 6 | AQUA CAPS MICRO SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 200g | |
| | Ref. 040CRMDR8 | • 6 | ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR Ø1,6 mm) - 1100g | |
| 6 | Ref. 042AMM | • 2 | MONO OT BOX MICRO SIZE (CASTABLE HOUSINGS FOR CAPS) | |
| | Ref. 041CAM | • 2 | STAINLESS STEEL HOUSINGS MICRO SIZE | |
| | Ref. 041CTM | • 2 | TITANIUM HOUSINGS MICRO SIZE | |
| | Ref. 101PPM | • 2 | PLASTIC POSITIONERS MICRO SIZE | |
| | Ref. 044PPM | • 2 | MICRO SPHERE ANALOGS | |
| | Ref. 100SAM | • 10 | ALUMINIUM SPACERS FOR OT BOX BARS MICRO SIZE | |

OT CAP ASSORTMENT SET



KIT NORMAL

Ref. 074KAN

KIT MICRO

Ref. 074KAM

SET CONTENT

- 8
- 8
- 6
- CONTENT

 CLEAR CAPS STANDARD RETENTION
 PINK CAPS SOFT RETENTION
 YELLOW CAPS EXTRA SOFT RETENTION
 GREEN CAPS VERY ELASTIC RETENTION
 GOLD CAPS (EXTRA RESILIENT) SLIGHTLY ELASTIC
 SILVER CAPS (EXTRA RESILIENT) ELASTIC AND GUMMY
 LIGHT GREEN CAPS STANDARD (FOR WORN OUT SPHERE)
 ORANGE CAPS ELASTIC (FOR WORN OUT SPHERE)
 BLACK CAPS LABORATORY USE ONLY
 TITAN CAPS

- TITAN CAPS



ASSORTMENT SET NORMAL & MICRO

SET CONTAINS

- 4+4 CLEAR CAPS STANDARD RETENTION
 4+4 PINK CAPS SOFT RETENTION
 2+2 YELLOW CAPS EXTRA SOFT RETENTION
 3+3 GREEN CAPS VERY ELASTIC RETENTION
 3+3 GOLD CAPS (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 3+3 SILVER CAPS (EXTRA RESILIENT) ELASTIC AND GUMMY

 2+2 LIGHT GREEN CAPS STANDARD (FOR WORN OUT SPHERE)

 2+2 ORANGE CAPS ELASTIC (FOR WORN OUT SPHERE)

 2+2 BLACK CAPS LABORATORY USE ONLY

- 1+1 TITAN CAPS

OT CAP ASSORTED RETENTIVE CAPS •

Ref. 074KNM



NORMAL

Ref. 064ACN


MICRO

Ref. 064ACM

SET CONTENT

- ON LENI
 GREEN CAPS VERY ELASTIC RETENTION
 YELLOW CAPS EXTRA SOFT RETENTION
 PINK CAPS SOFT RETENTION
 CLEAR CAPS STANDARD RETENTION
- 6
- 6

OT CAP SYSTEM ATTACHMENTS •



Ref. 038STF

TITANIUM FLEX SINGLE SPHERE NORMAL SIZE + TIN

FOR LASER SOLDERING - FOR BONDING

- SINGLE TITANIUM FLEX SPHERES PINK CAPS SOFT RETENTION SPHERE HOLDER

- BASI GUIDA CALCINABILI



Ref. 038STN

Ref. 038STM

TITANIUM SINGLE SPHERE NORMAL/MICRO SIZE + TIN

FOR LASER SOLDERING - FOR BONDING

SET CONTENT

- SINGLE TITANIUM SPHERES • 2
- PINK CAPS SOFT RETENTION SPHERE HOLDER
- CASTABLE SLIDING BASES

OT CAP NORMAL-MICRO TOOLS -

| | Ref. 485IC | 1 EXTRACTOR/INSERTER TOOL FOR CAPS OT EQUATOR - OT CAP NORMAL/MICRO |
|-----|-------------|--|
| | Ref. 74AC01 | 1 PARALLELOMETER MANDREL NORMAL/OT EQUATOR |
| K, | Ref. 74AC02 | 1 PARALLELOMETER MANDREL MICRO |
| | Ref. 74AC03 | 1 PARALLELOMETER MANDREL TECNO NORMAL AND MICRO SIZE (ONLY FOR OT CAP TECNO) |
| | Ref. 491EC | 1 MULTIUSE REMOVAL TOOL FOR RETENTIVE CAPS |
| | Ref. 008MDS | 1 SPHERE MEASURING GAUGE (TO CHECK THE SPHERICAL ATTACHMENT DIAMETERS: 1,6-1,8-2,2-2,5mm) |
| D. | Ref. 080RCN | • 1 REAMER TOOL FOR OT CAPS - NORMAL SIZE |
| A [| Ref. 080RCM | • 1 REAMER TOOL FOR OT CAPS - MICRO SIZE |
| | Ref. 082ATN | • 1 TOOL FOR TESTING CAP RETENTION - NORMAL SIZE |
| THE | Ref. 083ATM | • 1 TOOL FOR TESTING CAP RETENTION - MICRO SIZE |
| | | |

NORMAL - MICRO ACCESSORIES AND Ø2.2mm

| 1000° | Ref. 155CC | • 6 | CONNECTORS | |
|----------|----------------------------------|------|---|--|
| | Ref. 187CRS | • 2 | CONCAVE SPHERES AVAILABLE: Ø1.8 - 2.2 - 2.5mm | |
| | Ref. 187CRE | • 2 | CONCAVE OT EQUATOR | |
| ? | Ref. 189SRS | • 2 | SOLID SPHERES AVAILABLE: Ø1.8mm | |
| | Ref. 100PD | • 10 | PROTECTIVE DISKS (internal diameter 1,6mm) | |
| | Ref. 100PDR | • 10 | PINK PROTECTIVE DISKS (internal diameter 2,9mm) | |
| HEN | Ref. 161SCMN MICRO Ref. 161SCMM | • 4 | OT CAP NORMAL/MICRO SPACERS FOR HOUSING | |
| | Ref. 038GCS | • 2 | CASTABLE SLIDING BASES | |
| | Ref. 044PP22 | • 2 | MICRO ANALOGS Ø2.2mm | |

OT STRATEGY SPARE PARTS



Ref. 154PCS

4 OT STRATEGY STANDARD MALES



Ref. 054PLS

4 OT STRATEGY HIGH MALES



Ref. 154STY

2 CASTABLE STEADY



Ref. 047CSD AY

STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE (EXTRA SOFT RETENTION) - 500g



Ref. 047CSD SN

STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE (SOFT RETENTION) - 900g



Ref. 047CSD

STRATEGY CLEAR CAPS FOR DUPLICATION TECHNIQUE (STANDARD RETENTION) - 1300g



Ref. 047CLS

STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE (FOR LABORATORY USE ONLY)



Ref. 047CSDDR8

AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



Ref. 045CSM AY

 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING (EXTRA SOFT RETENTION) - 500g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM SN

STRATEGY PINK CAPS FOR S/STEEL HOUSING (SOFT RETENTION) - 900g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CSM

STRATEGY CLEAR CAPS FOR S/STEEL HOUSING (STANDARD RETENTION) -1300g (TO BE USED ONLY WITH HOUSINGS)



Ref. 045CLS

STRATEGY BLACK CAPS FOR S/STEEL HOUSING (FOR LABORATORY USE ONLY)



Ref. 045CSMDR8

AQUA STRATEGY CAPS FOR S/STEEL HOUSING (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE)



Ref. 046CSM

2 STRATEGY STAINLESS STEEL HOUSINGS



Ref. 1010TS

 2 STRATEGY PLASTIC POSITIONERS (FOR SOLDERING STAINLESS STEEL HOUSINGS)



OT STRATEGY ASSORTMENT CAP SET FOR DUPLICATION TECHNIQUE

SET CONTENT

- YELLOW CAPS EXTRA SOFT RETENTION PINK CAPS SOFT RETENTION CLEAR CAPS STANDARD RETENTION • 4

OT STRATEGY ASSORTMENT CAP SET FOR S/STEEL HOUSINGS

- YELLOW CAPS EXTRA SOFT RETENTION PINK CAPS SOFT RETENTION CLEAR CAPS STANDARD RETENTION Ref. 045ACS

Ref. 044PPS • 2 STRATEGY ANALOGS

• 1 INSERTER/EXTRACTOR TOOL FOR STRATEGY CAPS **Ref. 486ICS**

• 1 PARALLELOMETER MANDREL FOR STRATEGY CAPS Ref. 75AC04

• 1 REAMER TOOL FOR STRATEGY CAPS Ref. 081RCS

OT BAR MULTIUSE SPARE PARTS

Ref. 047ACS

• 2 CASTABLE BARS OT BAR MULTIUSE **Ref. 0220BB** 4 CASTABLE BOXES OT BAR MULTIUSE Ref. 025CPB • 4 POSITIONER CLIPS A Ref. 023CPA • 4 POSITIONER CLIPS B Ref. 024CPB

Ref. 027CRG

4 YELLOW CLIPS (MEDIUM RETENTION 1,8Kg)

Ref. 026CRR • 4 PINK CLIPS (SOFT RETENTION 1,2Kg)

Ref. 027CRS • 4 RED CLIPS (EXTRA SOFT RETENTION 0,6Kg)

Ref. 022RGO 3 CASTABLE CONNECTORS

• 1 INSERTER/EXTRACTOR TOOL FOR OT BAR MULTIUSE CLIPS **Ref. 429IOBM**

• 1 PARALLELOMETER KEY FOR OT BAR MULTIUSE **Ref. 0280CP**

OT VERTICAL SPARE PARTS -



OT REVERSE 3 SPARE PARTS

| | Ref. 120MR2 | •2 | RETENTIVE MALES IN TITANIUM + NYLON (COLOR CLEAR/ STANDARD RETENTION) | |
|---|-------------|-----|---|--|
| | Ref. 041CAM | • 2 | S/STEEL HOUSINGS MICRO SIZE | |
| | Ref. 121AAR | •2 | INOX ROOTS ANALOGS | |
| Ţ | Ref. 121AM | • 2 | PLASTIC MANUAL TOOLS | |
| • | Ref. 121PR | • 2 | PLASTIC POSITIONERS FOR DUPLICATION | |
| 9 | Ref. 100PDR | •10 | PINK PROTECTIVE DISKS | |
| | Ref. 124IMR | •1 | METAL INSERTING TOOL FOR OT REVERSE/3 MALE | |
| 0 | Ref. 72CBR | •1 | S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR | |
| 0 | Ref. 72CFR | • 1 | BRASS PARALLELOMETER KEY FORTITANIUM FEMALE | |
| | | | | |

OT LOCK SPARE PARTS



Ref. 882CS

SLIDING KEY (ONLY FOR TITANIUM - CASTABLE LOCKING PIN)



Ref. 882AF

• 1 FLEXIBLE RING



Ref. 882CB

• 1 LOCKING KEY



Ref. 882CPT

• 1 TITANIUM PREFABRICATED HOUSING (ONLY FOR TITANIUM - TITANIUM ADJUSTABLE LOCKING PIN)



Ref. 882PPC

• 1 CERAMIC PIN



Ref. 882CLSCT

• 1 ADJUSTABLE SLIDING KEY (ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)



Ref. 882DCC

CASTABLE SPACER RINGS (ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)

OT LOCK ACCESSORIES



Ref. 882CG

• 1 CONICAL GUIDE



Ref. 882CAS

• 1 UNLOCKING TOOL

SPHERE REPAIR SETS



Ø DISPONIBILI: 1,8 - 2,2 - 2,5 Ref. 087CRS

CONCAVE REPAIR SPHERE

SET CONTENT

· 2 CONCAVE SPHERES IN TITANIUM - TIN COATED

PINK CAPS SOFT RETENTION INSERTING TOOL

GAUGE AND STRIP HOLDER



Ref. 087CRE

CONCAVE REPAIR OT EQUATOR

SET CONTENT

CONCAVE OT EQUATOR IN TITANIUM - TIN COATED PINK CAPS SOFT RETENTION INSERTING TOOL

• 2 • 2

GAUGE AND STRIP HOLDER



Ref. 089SRS

SOLID REPAIR SPHERE

SET CONTENT

SOLID SPHERES 1.8mm IN TITANIUM - TIN COATED • 2 • 2

PINK CAPS SOFT RETENTION PROTECTIVE DISKS

KEY

DENTURE BOXES -



Ref. 0P6

SET CONTENT

• 6 ASSORTED BOXES, COLOURS: BLUE - LIGHT BLUE - WHITE



Ref. 100AS

• 30 SILICONE RINGS

SMALL PARALLELOMETER



Ref. 00PB


SMALL PARALLELOMETER DEVICE WITH MODEL HOLDER BASE (FOR COURSES, ETC.)

OT LEVER -

UNIVERSAL LEVER



OT LEVERFOR REMOVAL OF PROSTHESES AND ALIGNERS



OT CEM -

COMPOSITE CEMENT



Ref. OC

SET CONTENT

- TWO-COMPONENT SYRINGE 5ml
- •10 AUTOMIX TIPS




Ref. TOC

• 10 AUTOMIX TIPS

BROKEN SCREW EXTRACTOR KIT

BROKEN SCREW EXTRACTOR SETS ARE MANUFACTURED FOR MOST IMPLANT BRANDS





Ref. 680

SET CONTENT

- CLAW REAMER BUR REVERSE CUTTING BUR MANUAL EXTRACTOR CENTERING DEVICE

Ref. 680FS

REVERSE CUTTING BUR • 1

Ref. 680FA*

CLAW REAMER BUR

*Specify type of implant and diameter

THE BROKEN SCREW EXTRACTOR SET MAY HAVE SPACE PROBLEMS IN THE REAR SECTION, WHICH COULD AFFECT IT'S PERFORMANCE. IT'S IMPORTANT TO ANALYZE THE SITUATION BEFORE USING THE SET.



Live courses

The surgical-prosthetic process in all its phases, an excellent tool to facilitate learning new protocols.



Research and **Development**

Innovation and research is developed with the collaboration of numerous universities, both Italian and foreign.



Academy Rhein83



Conferences, events and practical courses, welco-ming professionals from all over the world for training in fixed and removable prosthetics.









DIRECT SERVICES

A telephone service is available during office hours, it furnishes direct answers for technical questions. You can address to Rhein83 for knowing the address of foreign distributors

Tel. +39 051 244510 - +39 051 244396 - Fax +39 051 245238

Rhein83 website www.rhein83.com

General informations nick@rhein83.it

Sales department international.orders@rhein83.it

Technical support direzione.tecnica@rhein83.it

Marketing department direzione@rhein83.it





Via ZAGO, 10/ABC 40128 - BOLOGNA (ITALY)

Tel. (+39) 051 244510 - (+39) 051 244396 Fax (+39) 051 245238







